


Visit our patient web page

[Learn more](#)

**1 Getting to know your mask**



F&P Solo™ Brochure



[Download](#)



Introduction to the F&P Solo

[Watch](#)

**2 Getting started**



**F&P myMask™ App**

With easy access to information and support, as you need it, F&P myMask App is the must-have mask companion. Access step-by-step videos of common tasks such as mask fitting, fine-tuning, cleaning and assembly.

**Scan to Download**



**3 Additional support material**



Fitting your mask

[Watch](#)



Cleaning your mask

[Watch](#)



User Instructions

[Download](#)

**Replacement schedule**

**Cushion**



Replace every 14 days or if the cushion begins to leak constantly, becomes discolored or is torn.

**Frame**



Replace every 3 months if the frame becomes discolored or cracked.

**Headgear**



Replace every 6 months or if the headgear becomes discolored, stretches beyond adjustment or the hook and fasteners no longer latch.