












SPECIFICATIONS

PRODUCT	 OJR410	 OJR412	 OJR414	 OJR416	 OJR418	 OJR520
NASAL CANNULA SIZE	XS	S	M	L	XL	XXL
F&P WigglewING™	 WJR210		 WJR212			 WJR214
F&P WigglewING						
Intended use	The F&P WigglewING™ is intended to secure a nasogastric (NG) tubing to the face in conjunction with an approved F&P respiratory interface					
Quantity	Box of 10					
Box components	F&P WigglewING, user instructions					
Weight	1.05 g		1.87 g			3.38 g
Disposal	Incineration or according to hospital protocol for adaptors					
PERFORMANCE						
Ambient range °C	18 – 26 °C / 64 – 79 °F					
Usage	Single patient use					
Duration of use	7 days					
Shelf life	3 years					
Compatible with	F&P Optiflow™ Junior 2 nasal cannula and F&P Optiflow™ Junior 2+ nasal cannula					
COMPONENTS AND COMPOSITIONS						
Predominant materials	Thermoplastic Elastomer; Hydrocolloid					
Materials not present	Not manufactured with natural rubber latex, PVC or Phthalates (DEHP, DBP, BBP)					
Manufacturing mode	Non-invasive device; Produced in a controlled working environment					
Disposal	Incineration or according to hospital protocol for medical waste					
REGULATORY						
Classification	AU-I, EU-IIa, Canada-I, US-I					
Country of origin	New Zealand					
Notified body	TÜV SÜD Product Service GmbH CE0123					

608131 REV B © 2019 Fisher & Paykel Healthcare Limited

Please note that the information in this specifications sheet (including product information and images) is summarized and provided for illustrative purposes only. Please refer to the relevant user instructions for more information and confirm details with your local Fisher & Paykel Healthcare representative prior to placing an order. Information subject to change without notice.