

NEW RT047 vented hospital full face mask with anti-asphyxiation valve



In response to our customers and patient's needs, Fisher & Paykel Healthcare is pleased to announce the release of a new vented version of our F&P Nivairo™ mask, expanding our range of full face masks.

Venting incorporated into the mask

The vented feature of the new Nivairo RT047 is designed to eliminate the need for an exhalation port. This helps to create an easier setup experience for healthcare professionals.


Venting designed to be quiet and comfortable

The unique venting configuration has been designed to reduce noise and direct draft away from the caregiver and patient's face.

Product order codes


Product code	Description	Quantity
RT047XS	Nivairo RT047 Extra Small	10/box
RT047S	Nivairo RT047 Small	10/box
RT047M	Nivairo RT047 Medium	10/box
RT047L	Nivairo RT047 Large	10/box

F&P Nivairo RT047 and Respironics AF541 venting leak rate*



- Respironics AF541
- F&P Nivairo RT047

* Leak rates of RT047 and AF541 without built-in exhalation port. AF541 leak rate data obtained through independent testing. RT047 data obtained through F&P internal testing.



Nivairo is a trademark of Fisher & Paykel Healthcare Limited. For patent information, see www.fphcare.com/ip

616845 REV A © 2019 Fisher & Paykel Healthcare Limited

Please contact your local representative for further information about this or any other Fisher & Paykel Healthcare product or therapy, or visit www.fphcare.com