

Breathing Circuits (Adult)

MODEL NUMBERS:

RT380: Adult Evaqua 2

RT385: Adult Evaqua 2 with pressure line

PERFORMANCE SPECS

Resistance to flow at 45L/min as packaged

Inspiratory	2.18 cm H ₂ O
Expiratory	1.83 cm H ₂ O
Length of tubing	1.5m
Compliance	2.1 mL/cm H ₂ O
Ambient Range	18-26C °C / 64-79 °F
Compressible Volume	1.60 L
Maximum chamber operating pressure	8 kPA
Duration of Use	14 days
Compatibility	Compatible with 700 series and MR850 humidifiers

FILTER

Filtration Efficiency

Viral	≥ 99.99%
Organism	ψ174 Bacteriophage
Bacterial	>99.999%
Organism	Bacillus subtilis
Mean particle size	3 µm

COMPONENTS & COMPOSITION

Filter included	Yes
Predominant materials	Polyethylene
Materials	Contains no natural rubber latex, may contain BPA or DEHP.
Manufacturing mode	Produced in a Controlled Working Environment
Disposal	Incineration or according to hospital protocol
Carton size	Box of 10



REGULATORY

Classification	Au Ila, EU Ila, Canada II, USA II
Country Of Origin	New Zealand
Intended Use	Single patient use, max 14 days
Notified Body	TÜV SÜD Product Services GmbH, 0123